

## CryoCool 8000 White

# Low temperature storage and liquid nitrogen labels

CryoCool 8000 White

### Thermal Transfer

#### Description:

CryoCool 8000 White is a 58 micron thermal transfer printable polypropylene label coated with a high performance acrylic adhesive for use in cryogenic applications. This adhesive offers resistance to temperatures as low as -196°C for liquid nitrogen applications

It is designed for use with Zebra 5095 resin thermal transfer ribbon to produce excellent print quality. The smear/scratch resistance and low temperature adhesion performance make CryoCool 8000 White an excellent choice for applications requiring a combination of image durability and extreme low temperature adhesion.

#### Suggested Applications:

Cryogenic applications involving a deep freezing process that takes objects down in temperature below -151 °C . CryoCool 8000 White will also withstand conditions such as dry ice (-80 C), steam autoclave and gamma radiation.

- Medical labs
- Universities/research facilities
- Hospitals
- Cold temperature/ industrial manufacturing

#### Technical Specifications

Description		Caliper
<b>Facestock</b>	White top coated polypropylene	58 microns
<b>Adhesive</b>	High Performance Permanent acrylic	20 microns
<b>Liner</b>	White Kraft liner	58 microns
<b>Total</b>		136 microns

**Recommended Zebra Ribbons:** 5095

**Minimum Application Temperature:** -29°C

**Service Temperature Range:** -196°C to 90°C

**Recommended Storage Conditions:** One year duration when stored at 0°C to 21°C



**CryoCool 8000 White**

(32°F to 70°F) at 35% to 50% RH

**Expected Exterior Life:** 6 months

**Test Procedure**

Labels were applied to glass vials (2.8 cm OD), polypropylene centrifuge tubes (3.5 cm OD, 50ml) and glass microscope slides and allowed a 24 hour dwell time before exposure to above conditions.

Environment	Test Method	Typical Results
High Temp.	30 days at listed temperature	No visible effect at 90 C (194 F)
Low Temp.	30 days at -70 C (-94 F)	No Visible effect
Freezer	3 cycles of 16 hours at -70 C (-94 F)/ 8 hours at room temp.	Glass vial : Recommended PP centrifuge tube: Recommended Glass microscope slide: Recommended Flat PP: Recommended
Pressure Cooker	3 cycles of 1 hour in 121 C (250 F) 15 psi pressure cooker/ 23 hours room temperature	Glass vial : Recommended PP centrifuge tube: Recommended Glass microscope slide: Recommended Flat PP: Recommended
Liquid Nitrogen	3 cycles of 4 hours at -196 C (-320 F)/ 20 hours at room temperature	Glass vial : Not recommended PP centrifuge tube: Recommended Glass microscope slide: Recommended Flat PP: Recommended
Freezer to boiling water	1 hour at -70 C (-94 F) then placed in boiling water 100 C (212 F)	Glass vial : May work, must test PP centrifuge tube: Recommended Glass microscope slide: May work, must test Flat PP: Recommended
Liquid Nitrogen to boiling water	1 hour at -196 C (-320 F) then placed in boiling water 100C (212 F) for 10 minutes	Glass vial : Not Recommended PP centrifuge tube: Recommended Glass microscope slide: May work, must test Flat PP: Recommended

**180 Degree Peel Room Temperature Peel Adhesion**

180 Peel (g/in)	Steel		Polyethylene		Polypropylene	
	5 min	24 hour	Loop Tack	24 hour	Loop Tack	24 hour
	810	908	360	408	317	408

**Product Performance and Suitability**

All information on this document is to be used for guidance only and is not to be used for setting specifications. All purchasers of Zebra products shall be responsible for independently determining if the product conforms to all requirements of the application.

